

You can put your trust in our system

- System is capable of handling up to 16 chairs
- Up to 3 years service life based on usage
- Trusted by dental colleges and military facilities
- Captures soluble and solid mercury per EPA test protocols
- Written guarantee of compliance
- Complimentary recycling **
- Reduce your mercury footprint
- No exposure to methylmercury
- Full warranty coverage during service life

** With installation of a new and/or replacement LibertyBOSS



Mars Bio-Med provides a UN/DOT approved storage and shipping container with prepaid return shipping and a certificate of recycling for the generator.

For the Recycling of:

- Chair-side traps
- Teeth with Amalgam
- Spent Amalgam Capsules
- Contact & Non-Contact Amalgam



CONTACT US

Toll-Free Canada & USA: 1.866.594.3648

Canada: Head Office, 1801 Wentworth St,
Unit 1, Whitby Ontario, L1N 8M2

USA: Warehouse, Depew, NY

info@marsbiomed.com
www.marsbiomed.com



MARS™

BIO-MED PROCESSES INC

The T-MACC System

The First Pre-Treatment System for Dental Office

Protecting our water from
Amalgam contamination

Toll-Free Canada & USA
1.866.594.3648
www.marsbiomed.com

Introducing the T-MACC System

Liberty BOSS

- Performance tested, independently verified - meets all EPA testing requirements
- No maintenance or filter changes; will save you money!
- Uses an ion exchange media that traps both solid and soluble mercury according to POTW and EPA test protocols
- Service life of up to 3 years based on practice actual use
- Largest sludge capacity on the market with 3,800 mL
- We use 1½" diameter internal lines to ensure clog-free operation



LINZ 3 Pack

- PowerLINZ daily cleaner and CleanLINZ quarterly cleaner
- Easy and safe to use
- Non-foaming, pH neutral
- Cost effective
- Anti-bacterial and Anti-microbial
- Environmentally friendly, 100% biodegradable
- Safe for all pump systems
- Compliant with CDC BMPs for cross contamination



The MARS Super Trap

The MARS Super Trap uses advanced technology as a pre-treatment solution to reduce your mercury footprint without modification to your existing equipment.

The MARS Super Trap replaces your existing chair side trap and works to reduce mercury and prophy buildup in the vacuum system.

Most popular trap numbers:

- #5500 pst
- #5501 pst
- #5503 pst
- #5505 pst
- #5512 pst

* additional sizes available

